


Anti-Rnf40 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O75150

Item Number: CM0005876



Introduction

Rabbit monoclonal antibody against human RNF40/BRE1B, validated for Western blot, IHC-P, and ELISA. Reacts with human and mouse. Useful for studying histone H2B monoubiquitination, epigenetic regulation, and KSHV lytic cycle.

[Learn More](#)

Parameter	Value
Product Name	RNF40 Rabbit mAb
Remarks/Alias	BRE1B; RBP95; STARING; RNF40
Species	Human
Gene ID	9810
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 231-431 of human RNF40 (O75150)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	O75150
Protein Weight	114kDa
Shipping	Ice bag