


Anti-Rnf40 Rabbit Polyclonal Antibody For Wb And Elisa - O75150

Item Number: CM0027508



Introduction

RNF40 Rabbit Polyclonal Antibody validated for WB and ELISA. Targets human E3 ubiquitin-protein ligase BRE1B (O75150), suitable for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	RNF40 Rabbit pAb
Remarks/Alias	BRE1B; RBP95; STARING; RNF40
Species	Human
Gene ID (Human)	9810
Gene ID	9810
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human RNF40 (NP_055586.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O75150
Protein Weight	114kDa
Shipping	Ice bag