


Anti-Rnf7 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Ubf6

Item Number: CM0004145



Introduction

Rabbit monoclonal antibody targeting human RNF7 (SAG/ROC2/RBX2/CKBBP1). Validated for WB, IHC-P, ELISA with cross-reactivity to mouse and rat. Suitable for ubiquitination and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	RNF7 Rabbit mAb
Remarks/Alias	SAG; ROC2; rbx2; CKBBP1; RNF7
Species	Human
Gene ID (Human)	9616
Gene ID	9616
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-113 of human RNF7 (NP_055060.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBF6
Protein Weight	13kDa
Shipping	Ice bag