


Anti-Rnf7 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Ubf6

Item Number: CM0030429



Introduction

Rabbit polyclonal antibody targeting RNF7 (also known as ROC2, SAG). Suitable for WB, IF/ICC, and ELISA applications with human samples. Recognizes a 13kDa protein (UniProt: Q9UBF6) involved in cullin-RING ligase complexes.

[Learn More](#)

Parameter	Value
Product Name	RNF7 Rabbit pAb
Remarks/Alias	SAG; ROC2; rbx2; CKBBP1; RNF7
Species	Human
GeneID	9616
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-113 of human RNF7 (NP_055060.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9UBF6
Protein Weight	13kDa
Shipping	Ice bag