


Anti-Rnf8 Polyclonal Antibody For Wb, Elisa - 076064

Item Number: CM0030017



Introduction

Rabbit polyclonal antibody targeting human RNF8, validated for western blotting and ELISA. Cross-reacts with mouse and rat. Suitable for DNA damage signaling and ubiquitination research.

[Learn More](#)

Parameter	Value
Product Name	RNF8 Rabbit pAb
Remarks/Alias	hRNF8; RNF8
Species	Human
GeneID (human)	9025
GeneID	9025
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 230-485 of human RNF8 (NP_003949.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O76064
Protein Weight	56kDa
Shipping	Ice bag