


Anti-Rnf8 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O76064

Item Number: CM0006724



Introduction

RNF8 rabbit monoclonal antibody validated for Western blot, IHC-P, and ELISA. Detects human, mouse, and rat RNF8 E3 ubiquitin ligase (~56 kDa) involved in DNA damage repair.

[Learn More](#)

Parameter	Value
Product Name	RNF8 Rabbit mAb
Remarks/Alias	hRNF8; RNF8
Species	Human
GeneID (Human)	9025
GeneID	9025
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human RNF8 (NP_003949.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O76064
Protein Weight	56kDa
Shipping	Ice bag