


## Anti-Rngtt Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O60942

Item Number: CM0028820



### Introduction

Anti-RNGTT rabbit polyclonal antibody for western blot, IF/ICC, and ELISA. Reacts with human, mouse, and rat RNGTT. Ideal for research on mRNA cap formation and diagnostic assay development.

[Learn More](#)

| Parameter        | Value                                                                                                                                 |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | RNGTT Rabbit pAb                                                                                                                      |
| Remarks/Alias    | HCE; HCE1; hCAP; CAP1A; RNGTT                                                                                                         |
| Species          | Human                                                                                                                                 |
| GeneID (Human)   | 8732                                                                                                                                  |
| GeneID           | 8732                                                                                                                                  |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human RNGTT (NP_003791.3). |
| Source           | Rabbit                                                                                                                                |
| Category         | Polyclonal Antibodies                                                                                                                 |
| Application      | WB, IF/ICC, ELISA                                                                                                                     |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                     |
| SWISS            | O60942                                                                                                                                |
| Protein Weight   | 69kDa                                                                                                                                 |
| Shipping         | Ice bag                                                                                                                               |