

Anti-Rnh1 Polyclonal Antibody For Wb, If/Icc, Elisa - P13489

Item Number: CM0026172



Introduction

High-quality rabbit polyclonal antibody against RNH1 (Ribonuclease inhibitor, RAI). Validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	RNH1 Rabbit pAb
Remarks/Alias	RAI; RNH; RNH1
Species	Human
GeneID (Human)	6050
GeneID	6050
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 202-461 of human RNH1 (NP_002930.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13489
Protein Weight	50kDa
Shipping	Ice bag