

Anti-Robo1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y6N7

Item Number: CM0018535



Introduction

Rabbit polyclonal antibody against human ROBO1, suitable for WB, IF/ICC, and ELISA applications. Cross-reacts with mouse and rat. Ideal for axon guidance, neurodevelopment, and cell migration studies.

[Learn More](#)

Parameter	Value
Product Name	ROBO1 Rabbit pAb
Remarks/Alias	NORS; NYS8; SAX3; CPHD8; DUTT1; ROBO1
Species	Human
GeneID (human)	6091
GeneID	6091
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1452-1651 of human ROBO1 (NP_002932.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6N7
Protein Weight	181kDa
Shipping	Ice bag