

Anti-Rock1 Rabbit Monoclonal Antibody For Wb, Elisa - Q13464

Item Number: CM0002699



Introduction

ROCK1 Rabbit Monoclonal Antibody validates in WB and ELISA with reactivity across human, mouse, and rat. Ideal for cytoskeleton dynamics, cell migration, and signaling research applications.

[Learn More](#)

Parameter	Value
Product Name	ROCK1 Rabbit mAb
Remarks/Alias	ROCK-I; P160ROCK; ROCK1
Species	Human
Gene ID	6093
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1100-1200 of human ROCK1 (Q13464).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13464
Protein Weight	158kDa
Shipping	Ice bag