

Anti-Rock1 Rabbit Polyclonal Antibody For Wb, Ip, If, Elisa - Q13464

Item Number: CM0014899



Introduction

ROCK1 rabbit polyclonal antibody validated for WB, IP, IF-P, ELISA in human, mouse, rat samples. Key regulator of actin cytoskeleton and cell polarity, ideal for cell motility and cancer research.

[Learn More](#)

Parameter	Value
Product Name	ROCK1 Rabbit pAb
Remarks/Alias	ROCK-I; P160ROCK; ROCK1
Species	Human
GeneID (Human)	6093
GeneID	6093
Immunogen	Recombinant protein corresponding to amino acids 400-660 of human ROCK1 (NP_005397.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13464
Protein Weight	158kDa
Shipping	Ice bag