

## Anti-Rock1 Rabbit Polyclonal Antibody (Ko Validated) For Wb, Ihc-P, Elisa - Q13464

Item Number: CM0015011



### Introduction

High-quality rabbit polyclonal antibody against human ROCK1, validated by knockout, for Western blot, IHC-P, and ELISA. Reacts with human and mouse samples. Ideal for cytoskeleton, cell polarity, and signaling research.

[Learn More](#)

| Parameter        | Value                                                                                                                                  |
|------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | [KO Validated] ROCK1 Rabbit pAb                                                                                                        |
| Remarks/Alias    | ROCK-I; P160ROCK; K1                                                                                                                   |
| Species          | Human                                                                                                                                  |
| GeneID (Human)   | 6093                                                                                                                                   |
| GeneID           | 6093                                                                                                                                   |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 400-660 of human ROCK1 (NP_005397.1) |
| Source           | Rabbit                                                                                                                                 |
| Category         | Polyclonal Antibodies                                                                                                                  |
| Application      | WB, IHC-P, ELISA                                                                                                                       |
| Cross Reactivity | Human, Mouse                                                                                                                           |
| SWISS            | Q13464                                                                                                                                 |
| Protein Weight   | 158kDa                                                                                                                                 |
| Shipping         | Ice bag                                                                                                                                |