

Anti-Rock2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O75116

Item Number: CM0027442



Introduction

Rabbit polyclonal antibody against ROCK2, validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat ROCK2 at ~161kDa. Ideal for cytoskeleton and kinase signaling research.

[Learn More](#)

Parameter	Value
Product Name	ROCK2 Rabbit pAb
Remarks/Alias	ROCK-II; ROCK2
Species	Human
Gene ID (Human)	9475
Gene ID	9475
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1109-1388 of human ROCK2 (NP_004841.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75116
Protein Weight	161kDa
Shipping	Ice bag