


Anti-Rogdi Polyclonal Antibody For Wb, Elisa - Q9Gzn7

Item Number: CM0025354



Introduction

Rabbit polyclonal antibody targeting human ROGDI, validated for Western blot (WB) and ELISA, with cross-reactivity to mouse and rat. Suitable for neuroscience research and disease model studies.

[Learn More](#)

Parameter	Value
Product Name	ROGDI Rabbit pAb
Remarks/Alias	KTZS; RAV2; ROGDI1; ROGDI
Species	Human
Gene ID (Human)	79641
Gene ID	79641
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human ROGDI (NP_078865.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9GZN7
Protein Weight	32kDa
Shipping	Ice bag