

Anti-Romo1 Rabbit Monoclonal Antibody For Wb - P60602

Item Number: CM0010738

Introduction

Rabbit monoclonal antibody specifically targeting human ROMO1, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for ROS and mitochondrial morphology studies.

[Learn More](#)

Parameter	Description
Product Name	[KD Validated] ROMO1 Rabbit mAb
Remarks/Alias	MTGM, MTGMP, C20orf52, bA353C18.2, ROMO1
Species	Human
GeneID (Human)	140823
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-79 of human ROMO1 (NP_542786.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-F, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60602
Protein Weight	8kDa
Shipping	Ice bag