


Anti-Ror1 Polyclonal Antibody For Wb, Elisa - Q01973

Item Number: CM0028354



Introduction

ROR1 Rabbit Polyclonal Antibody validated for Western blot and ELISA, cross-reactive with rat. Ideal for Wnt signaling, neurotrophic receptor studies, and cancer research. Unconjugated, 104kDa target.

[Learn More](#)

Parameter	Value
Product Name	ROR1 Rabbit pAb
Remarks/Alias	NTRKR1; dJ537F10.1; ROR1
Species	Human
Gene ID (Human)	4919
Gene ID	4919
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-403 of human ROR1 (NP_005003.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q01973
Protein Weight	104kDa
Shipping	Ice bag