


## Anti-Ror2 Rabbit Monoclonal Antibody For Wb And Elisa - Q01974

Item Number: CM0012048



### Introduction

Rabbit monoclonal antibody targeting human ROR2, suitable for Western blot and ELISA. Cross-reacts with mouse. Useful for studies on chondrogenesis, osteogenesis, and Wnt signaling. Detects ~105 kDa protein.

[Learn More](#)

| Parameter        | Value                                                                                            |
|------------------|--------------------------------------------------------------------------------------------------|
| Product Name     | ROR2 Rabbit mAb                                                                                  |
| Remarks/Alias    | BDB; BDB1; NTRKR2; ROR2                                                                          |
| Species          | Human                                                                                            |
| Gene ID (Human)  | 4920                                                                                             |
| Gene ID          | 4920                                                                                             |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ROR2 (Q01974). |
| Source           | Rabbit                                                                                           |
| Category         | Monoclonal Antibodies                                                                            |
| Application      | WB, ELISA                                                                                        |
| Cross Reactivity | Mouse                                                                                            |
| SWISS            | Q01974                                                                                           |
| Protein Weight   | 105kDa                                                                                           |
| Shipping         | Ice bag                                                                                          |