

Anti-Ror2 Rabbit Polyclonal Antibody For Wb, Elisa - Q01974

Item Number: CM0024771



Introduction

Rabbit polyclonal antibody targeting human ROR2, validated for WB and ELISA. Recognizes human, mouse, and rat ROR2. Ideal for research on Wnt signaling, osteogenesis, and cartilage development.

[Learn More](#)

Parameter	Value
Product Name	ROR2 Rabbit pAb
Remarks/Alias	BDB; BDB1; NTRKR2; ROR2
Species	Human
GeneID	4920
Immunogen	Synthetic peptide (aa 100-200 of human ROR2, NP_004551.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01974
Protein Weight	105kDa
Shipping	Ice bag