


Anti-Rora Polyclonal Antibody For Wb, Elisa - P35398

Item Number: CM0026291



Introduction

RORA rabbit polyclonal antibody (P35398) validated for Western blot and ELISA. Detects human, mouse, and rat RORA. Ideal for circadian rhythm, metabolism, and T-cell differentiation research.

[Learn More](#)

Parameter	Value
Product Name	RORA Rabbit pAb
Remarks/Alias	ROR1; ROR2; ROR3; RZRA; NR1F1; RORa1; IDDECA; RORalpha; RZR-ALPHA; RORA
Species	Human
GeneID (Human)	6095
GeneID	6095
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 175-255 of human RORA (NP_599023.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35398
Protein Weight	53KDa/59kDa/61KDa/63KDa
Shipping	Ice bag