

Anti-Rorc Polyclonal Antibody For Wb, Elisa - P51449

Item Number: CM0013323



Introduction

Rabbit polyclonal anti-RORC antibody for WB and ELISA. Detects human RORC, cross-reacts with mouse and rat. Ideal for studying nuclear receptor ROR-gamma in immune cell differentiation, circadian rhythm, and metabolic regulation.

[Learn More](#)

Parameter	Value
Product Name	RORC Rabbit pAb
Remarks/Alias	TOR; RORG; RZRG; IMD42; NR1F3; RZR-GAMMA; RORC
Species	Human
Gene ID (Human)	6097
Gene ID	6097
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 223-312 of human RORC (NP_005051.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51449
Protein Weight	58kDa
Shipping	Ice bag