

## Anti-Rorβ Rabbit Monoclonal Antibody For Wb, Elisa - Q92753

Item Number: CM0002512

### Introduction

Validated rabbit monoclonal antibody against human RORβ for Western blot and ELISA. Reacts with human, mouse, rat. Targets nuclear receptor ROR-beta, essential for circadian rhythm, vision, and neural development.

[Learn More](#)

| Parameter        | Value                                                   |
|------------------|---------------------------------------------------------|
| Product Name     | RORβ Rabbit mAb                                         |
| Remarks/Alias    | RZRB; EIG15; NR1F2; ROR-BETA; RZR-BETA; bA133M9.1; RORβ |
| Species          | Human                                                   |
| GeneID (Human)   | 6096                                                    |
| GeneID           | 6096                                                    |
| Immunogen        | Synthetic Peptide                                       |
| Source           | Rabbit                                                  |
| Category         | Monoclonal Antibodies                                   |
| Application      | WB, ELISA                                               |
| Cross Reactivity | Human, Mouse, Rat                                       |
| SWISS            | Q92753                                                  |
| Protein Weight   | 53kDa                                                   |
| Shipping         | Ice bag                                                 |