


Anti-Rpa2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P15927

Item Number: CM0022386



Introduction

Rabbit polyclonal antibody against human RPA2 (RP-A p32), validated for WB, IHC-P, IF/ICC, ELISA. Essential for DNA replication and DNA damage response. Ideal for studying DNA repair and genome stability.

[Learn More](#)

Parameter	Value
Product Name	RPA2 Rabbit pAb
Remarks/Alias	RPA2; REPA2; RP-A p32; RP-A p34; RPA32; RP-Ap32; RP-Ap34; replication protein A2
Species	Human
GeneID (Human)	6118
GeneID	6118
Immunogen	Recombinant protein corresponding to a sequence within amino acids 1-270 of human RPA2(NP_002937.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P15927
Protein Weight	29kDa/30kDa/38kDa
Shipping	Ice bag