


Anti-Rpa32/Rpa2 Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - P15927

Item Number: CM0010715



Introduction

Rabbit monoclonal RPA32/RPA2 antibody validated for WB, IF/ICC, IHC-P, ELISA in human, mouse, rat. Essential for DNA replication and damage repair. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	RPA32/RPA2 Rabbit mAb
Remarks/Alias	REPA2; RPA32; RP-A p32; RP-A p34
Species	Human
GeneID (Human)	6118
Immunogen	Recombinant protein corresponding to amino acids 1-270 of human RPA32/RPA2 (NP_002937.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15927
Protein Weight	29kDa
Shipping	Ice bag