

Anti-Rpa70/Rpa1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P27694

Item Number: CM0018051



Introduction

Rabbit polyclonal antibody against human RPA70/RPA1 (Replication protein A 70 kDa subunit), validated for Western blot, immunohistochemistry, and ELISA. Useful for studying DNA replication, repair, and damage response.

[Learn More](#)

Parameter	Value
Product Name	RPA70/RPA1 Rabbit pAb
Remarks/Alias	HSSB; RF-A; RP-A; REPA1; RPA70; MST075; PFBMFT6; RPA70/RPA1
Species	Human
GeneID (Human)	6117
GeneID	6117
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human RPA70/RPA1 (NP_002936.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P27694
Protein Weight	68kDa
Shipping	Ice bag