

Anti-Rpap2 Polyclonal Antibody For Wb And Elisa - Q8Ixx5

Item Number: CM0020294



Introduction

Rabbit polyclonal antibody targeting human RPAP2 (RNA polymerase II-associated protein 2). Suitable for WB and ELISA applications. RPAP2 is a CTD phosphatase involved in snRNA transcription and ER stress signaling. SWISS: Q8IXW5.

[Learn More](#)

Parameter	Value
Product Name	RPAP2 Rabbit pAb
Remarks/Alias	Rtr1; C1orf82; RPAP2
Species	Human
GeneID (Human)	79871
GeneID	79871
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-612 of human RPAP2 (Q8IXW5)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8IXW5
Protein Weight	70kDa
Shipping	Ice bag