


Anti-Rpb1 Ctd Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Chip, Elisa - P24928

Item Number: CM0008479



Introduction

High-quality rabbit monoclonal antibody targeting Rpb1 CTD, validated for WB, IF/ICC, IHC-P, ChIP, and ELISA. Detects human, mouse, and rat Rpb1 (~217 kDa). Ideal for studying RNA polymerase II transcription regulation. One-stop IVD sourcing from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	Rpb1 CTD Rabbit mAb
Remarks/Alias	RPB1; RPO2; POLR2; POLRA; RPBh1; RPOL2; NEDHIB; RpiILS; hsRPB1; hRPB220
Species	Human
Gene ID	5430
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1611-1710 of human Rpb1 CTD (NP_000928.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ChIP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24928
Protein Weight	217 kDa/64 kDa
Shipping	Ice bag