

Anti-Rpe Polyclonal Antibody For Wb - Q96At9

Item Number: CM0016228



Introduction

High-quality rabbit polyclonal antibody targeting human Ribulose-phosphate 3-epimerase (RPE, Q96AT9). Suitable for Western blotting, immunofluorescence, and ELISA. Shows reactivity with human, mouse, and rat samples. Detects an expected band at 25kDa.

[Learn More](#)

Parameter	Value
Product Name	RPE Rabbit pAb
Remarks/Alias	RPE2-1; RPE
Species	Human
GeneID (human)	6120
GeneID	6120
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-228 of human RPE (NP_954699.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96AT9
Protein Weight	25kDa
Shipping	Ice bag