

Anti-Rpe65 Polyclonal Antibody For Wb And Elisa - Q16518

Item Number: CM0022579

Introduction

Rabbit polyclonal antibody targeting human RPE65, validated for WB and ELISA. Cross-reacts with mouse. RPE65 is a critical isomerohydrolase in the retinoid cycle, essential for vision. Ideal for retinal disease research.

[Learn More](#)

Parameter	Value
Product Name	RPE65 Rabbit pAb
Remarks/Alias	p63; BCO3; LCA2; RP20; rd12; mRPE65; sRPE65; RPE65
Species	Human
GeneID (human)	6121
GeneID	6121
Immunogen	Recombinant protein corresponding to a sequence within amino acids 340-440 of human RPE65 (NP_000320.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q16518
Protein Weight	61kDa
Shipping	Ice bag