

Anti-Rpe65 Polyclonal Antibody For Wb, If-P, Elisa - Q16518

Item Number: CM0025690

Introduction

Rabbit anti-RPE65 polyclonal antibody validated for Western blot (WB), immunofluorescence paraffin (IF-P), and ELISA. Detects human, mouse, and rat retinoid isomerohydrolase; ideal for visual cycle and inherited retinal disease research.

[Learn More](#)

Parameter	Value
Product Name	RPE65 Rabbit pAb
Remarks/Alias	p63; BCO3; LCA2; RP20; rd12; mRPE65; sRPE65; RPE65
Species	Human
Gene ID	6121
Immunogen	Recombinant fusion protein corresponding to amino acids 274-533 of human RPE65 (NP_000320.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16518
Protein Weight	61kDa
Shipping	Ice bag