

Anti-Rpe65 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q16518

Item Number: CM0005333

Introduction

High-specificity rabbit monoclonal antibody against RPE65, validated for WB, IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for retinoid cycle research and vision studies.

[Learn More](#)

Parameter	Value
Product Name	RPE65 Rabbit mAb
Remarks/Alias	p63; BCO3; LCA2; RP20; rd12; mRPE65; sRPE65; RPE65
Species	Human
GeneID (Human)	6121
GeneID	6121
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 400-500 of human RPE65 (Q16518).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q16518
Protein Weight	61kDa
Shipping	Ice bag