


## Anti-Rpf2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H7B2

Item Number: CM0023401



### Introduction

RPF2 Rabbit Polyclonal Antibody validated for Western blot and ELISA. Detects human, mouse, rat RPF2, a nucleolar protein involved in ribosome large subunit assembly. Suitable for ribosome biogenesis research.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | RPF2 Rabbit pAb                                                                                                  |
| Remarks/Alias    | BXDC1; bA397G5.4; RPF2                                                                                           |
| Species          | Human                                                                                                            |
| GeneID (Human)   | 84154                                                                                                            |
| GeneID           | 84154                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-306 of human RPF2 (NP_115570.1). |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | Q9H7B2                                                                                                           |
| Protein Weight   | 36kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |