


Anti-Rpgr Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q92834

Item Number: CM0024606



Introduction

Rabbit polyclonal antibody targeting human RPGR (Q92834) for WB, IF-P, IHC-P, ELISA. Detects human, mouse, and rat RPGR, a GEF critical for ciliary trafficking and retinal homeostasis. Ideal for vision and ciliopathy research.

[Learn More](#)

Parameter	Value
Product Name	RPGR Rabbit pAb
Remarks/Alias	CRD; RP3; COD1; PCDX; RP15; XLRP3; orf15; CORDX1; RPGR
Species	Human
Gene ID (Human)	6103
Gene ID	6103
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human RPGR (NP_001030025.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92834
Protein Weight	113kDa
Shipping	Ice bag