

Anti-Rpgr Rabbit Monoclonal Antibody For Wb, If, Elisa - Q92834-6

Item Number: CM0006962

Introduction

High-quality rabbit monoclonal antibody targeting human RPGR (isoform 6). Validated for WB, IF/ICC, IF-P, and ELISA, with cross-reactivity to mouse and rat. Ideal for ciliary and retinal disease research.

[Learn More](#)

Parameter	Specification
Product Name	RPGR Rabbit mAb
Remarks/Alias	CRD; RP3; COD1; PCDX; RP15; XLRP3; orf15; CORDX1; RPGR
Species	Human
GeneID (Human)	6103
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human RPGR. (NP_001030025.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92834-6
Protein Weight	52kDa/70kDa/83kDa/90kDa/113kDa/127kDa
Shipping	Ice bag