

Anti-Rph3A (Rabphilin-3A) Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y2J0

Item Number: CM0029120



Introduction

Rabbit polyclonal antibody targeting human Rabphilin-3A (RPH3A, Exophilin-1). Validated for WB and ELISA; cross-reacts with mouse and rat. For studying regulated exocytosis and synaptic function.

[Learn More](#)

Parameter	Value
Product Name	RPH3A Rabbit pAb
Remarks/Alias	RPH3A
Species	Human
Gene ID	22895 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human RPH3A (NP_055769.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9Y2J0
Protein Weight	77kDa
Shipping	Ice bag