

Anti-Rph3Al Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Une2

Item Number: CM0026466

Introduction

Rabbit polyclonal antibody targeting RPH3AL (NOC2), validated for WB, IF/ICC, and ELISA. Detects human and mouse RPH3AL, a Rab effector protein crucial for regulated exocytosis in endocrine and exocrine cells.

[Learn More](#)

Parameter	Value
Product Name	RPH3AL Rabbit pAb
Remarks/Alias	NOC2; RPH3AL
Species	Human
Gene ID (Human)	9501
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 161-315 of human RPH3AL (NP_008918.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UNE2
Protein Weight	34kDa
Shipping	Ice bag