

Anti-Rpia Polyclonal Antibody For Wb, Ihc-P, Elisa - P49247

Item Number: CM0026717

Introduction

Rabbit polyclonal antibody targeting human RPIA (ribose-5-phosphate isomerase), validated for WB, IHC-P, ELISA with cross-reactivity in mouse and rat. Ideal for pentose phosphate pathway studies.

[Learn More](#)

Parameter	Value
Product Name	RPIA Rabbit pAb
Remarks/Alias	RPI; RPIAD; RPIA
Species	Human
Gene ID (Human)	22934
Gene ID	22934
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 82-311 of human RPIA (NP_653164.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49247
Protein Weight	33kDa
Shipping	Ice bag