

Anti-Rpl10A Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P62906

Item Number: CM0027989



Introduction

Rabbit polyclonal antibody targeting human RPL10A (ribosomal protein L10a), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for translation and ribosome studies.

[Learn More](#)

Parameter	Value
Product Name	RPL10A Rabbit pAb
Remarks/Alias	uL1; L10A; CSA19; NEDD6; Csa-19; RPL10A
Species	Human
Gene ID (Human)	4736
Gene ID	4736
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human RPL10A (NP_009035.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62906
Protein Weight	25kDa
Shipping	Ice bag