


Anti-Rpl11 Polyclonal Antibody For If/Icc, Elisa - P62913

Item Number: CM0015233



Introduction

CamelBio's anti-RPL11 rabbit polyclonal antibody targets human 60S ribosomal protein L11 (uL5). Validated for IF/ICC and ELISA, it is ideal for studying ribosome biogenesis, translational control, and p53 signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	RPL11 Rabbit pAb
Remarks/Alias	L11; uL5; DBA7; GIG34; RPL11
Species	Human
Gene ID (Human)	6135
Gene ID	6135
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-177 of human RPL11 (NP_001186731.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P62913
Protein Weight	20kDa
Shipping	Ice bag