


Anti-Rpl11 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P62913

Item Number: CM0005920



Introduction

Rabbit monoclonal antibody against human RPL11, suitable for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat species. Ideal for ribosome research and p53 pathway studies.

[Learn More](#)

Parameter	Value
Product Name	RPL11 Rabbit mAb
Remarks/Alias	L11, uL5, DBA7, GIG34
Species	Human
GeneID (Human)	6135
GeneID	6135
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-178 of RPL11 (NP_000966.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62913
Protein Weight	20kDa
Shipping	Ice bag