

Anti-Rpl11 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P62913

Item Number: CM0030103



Introduction

Rabbit polyclonal antibody targeting human RPL11 (large ribosomal subunit protein uL5), validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat. Suitable for ribosome biogenesis, protein synthesis, and p53 signaling studies.

[Learn More](#)

Parameter	Value
Product Name	RPL11 Rabbit pAb
Remarks/Alias	L11; uL5; DBA7; GIG34; RPL11
Species	Human
GeneID (Human)	6135
GeneID	6135
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-177 of RPL11 (NP_001186731.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62913
Protein Weight	20kDa
Shipping	Ice bag