

Anti-Rpl13 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P26373

Item Number: CM0005171



Introduction

Rabbit monoclonal antibody targeting human RPL13 (eL13) for Western blotting, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for ribosomal protein and bone development research.

[Learn More](#)

Parameter	Value
Product Name	RPL13 Rabbit mAb
Remarks/Alias	L13; BBC1; eL13; D16S44E; SEMDIST; D16S444E; RPL13
Species	Human
Gene ID (Human)	6137
Gene ID	6137
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human RPL13 (NP_150254.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26373
Protein Weight	18kDa/24kDa
Shipping	Ice bag