

Anti-Rpl13 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P26373

Item Number: CM0026837



Introduction

Rabbit polyclonal antibody targeting human RPL13 (60S ribosomal protein L13), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for ribosome and protein synthesis studies.

[Learn More](#)

Parameter	Value
Product Name	RPL13 Rabbit pAb
Remarks/Alias	L13; BBC1; eL13; D16S44E; SEMDIST; D16S444E; RPL13
Species	Human
GeneID (Human)	6137
GeneID	6137
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-211 of human RPL13 (NP_150254.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26373
Protein Weight	24kDa
Shipping	Ice bag