


Anti-Rpl13A Rabbit Polyclonal Antibody For Wb, Elisa - P40429

Item Number: CM0018767



Introduction

High-quality rabbit polyclonal antibody against human RPL13A (uL13). Validated for Western blot (WB) and ELISA. Reacts with human, mouse, rat. Ideal for ribosomal function, translation regulation, and GAIT complex studies.

[Learn More](#)

Parameter	Value
Product Name	RPL13A Rabbit pAb
Remarks/Alias	L13A; uL13; TSTA1; RPL13A
Species	Human
GeneID (Human)	23521
GeneID	23521
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-197 of human RPL13A (NP_036555.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40429
Protein Weight	24kDa
Shipping	Ice bag