


Anti-Rpl17 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P18621

Item Number: CM0026468



Introduction

Rabbit polyclonal antibody against human RPL17 (60S ribosomal protein L17, alias L17, PD-1, uL22) for WB, IHC-P, IF/ICC, ELISA. Validated in human, mouse, rat. Suitable for ribosomal protein research and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	RPL17 Rabbit pAb
Remarks/Alias	L17; PD-1; uL22; RPL23; RPL17
Species	Human
GeneID (Human)	6139
GeneID	6139
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-184 of human RPL17 (NP_000976.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18621
Protein Weight	21kDa
Shipping	Ice bag