


Anti-Rpl18 Polyclonal Antibody For Wb, Elisa - Q07020

Item Number: CM0013068



Introduction

Rabbit polyclonal antibody targeting human RPL18 (60S ribosomal protein L18). Validated for Western blot and ELISA, cross-reacts with mouse and rat. Suitable for studying ribosome biogenesis, protein synthesis, and Diamond-Blackfan anemia.

[Learn More](#)

Parameter	Value
Product Name	RPL18 Rabbit pAb
Remarks/Alias	L18; eL18; DBA18; RPL18
Species	Human
Gene ID (Human)	6141
Gene ID	6141
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-188 of human RPL18 (NP_000970.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07020
Protein Weight	22kDa
Shipping	Ice bag