


Anti-Rpl18A Polyclonal Antibody For Wb And Elisa - Q02543

Item Number: CM0023581



Introduction

Rabbit polyclonal antibody against human RPL18A (60S ribosomal protein L18a/eL20), validated for Western blot and ELISA. Detects endogenous levels in human, mouse, rat samples. Ideal for ribosomal biology and translational control research.

[Learn More](#)

Parameter	Value
Product Name	RPL18A Rabbit pAb
Remarks/Alias	L18A; eL20; RPL18A
Species	Human
Gene ID (Human)	6142
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-176 of human RPL18A (NP_000971.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02543
Protein Weight	21kDa
Shipping	Ice bag