

Anti-Rpl19 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P84098

Item Number: CM0023687



Introduction

Rabbit polyclonal antibody against human RPL19 (60S ribosomal protein L19). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for ribosomal and translation studies.

[Learn More](#)

Parameter	Value
Product Name	RPL19 Rabbit pAb
Remarks/Alias	L19; eL19; RPL19
Species	Human
GeneID (Human)	6143
GeneID	6143
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human RPL19 (NP_000972.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P84098
Protein Weight	23kDa
Shipping	Ice bag