


Anti-Rpl21 Rabbit Polyclonal Antibody - P46778

Item Number: CM0017868



Introduction

RPL21 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Recognizes human, mouse, rat RPL21. Useful for studying ribosomal function and protein synthesis.

[Learn More](#)

Parameter	Value
Product Name	RPL21 Rabbit pAb
Remarks/Alias	L21; eL21; HYPT12; RPL21
Species	Human
GeneID (human)	6144
GeneID	6144
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human RPL21 (NP_000973.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46778
Protein Weight	19kDa
Shipping	Ice bag