


Anti-Rpl23 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P62829

Item Number: CM0028856



Introduction

High-quality rabbit polyclonal antibody against human RPL23 (60S ribosomal protein L23), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. For ribosome biology studies.

[Learn More](#)

Parameter	Value
Product Name	RPL23 Rabbit pAb
Remarks/Alias	L23; uL14; rpl17; RPL23
Species	Human
GeneID (Human)	9349
GeneID	9349
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human RPL23 (NP_000969.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62829
Protein Weight	15kDa
Shipping	Ice bag