


Anti-Rpl23A Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P62750

Item Number: CM0026133



Introduction

Rabbit polyclonal antibody against human RPL23A (60S ribosomal protein L23a), validated for WB, IHC-P, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for ribosome research and cancer studies.

[Learn More](#)

Parameter	Specification
Product Name	RPL23A Rabbit pAb
Remarks/Alias	L23A; uL23; MDA20; RPL23A
Species	Human
GeneID (Human)	6147
GeneID	6147
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-156 of human RPL23A (NP_000975.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62750
Protein Weight	18kDa
Shipping	Ice bag